





















SatShare™ Model	Validated SatShare™ Monitor Interfaces	Model / Plug-In Parameter Modules	Host Instrument Connector
AT01 P/N 1595	Airshields Athena Monitor with Nellcor SpO <sub>2</sub>		
DO01 P/N 1324	Datex-Ohmeda AS/3 Patient Monitor	M-NSAT	
	System 5 Compact Monitor	M-NSAT	
DO02 P/N 1326	Datex-Ohmeda AS/3 Patient Monitor	M-NESTPR-00-00 M-NESTR-00-00	
	Datex-Ohmeda CardioCap/5 Monitor		
	Datex-Ohmeda Light Monitor		
	Datex-Ohmeda System 5 Compact Monitor	M-NESTPR-00-00 M-NESTR-00-00	
DO03 P/N 1528	Datex-Ohmeda AS/3 Patient Monitor	M-EST	
DO04 P/N 1533	Datex-Ohmeda AS/3 Patient Monitor	M-ESTP	
	Datex-Ohmeda System 5 Compact Monitor	M-EST M-ESTP	
DO05 P/N 2023	Datex-Ohmeda System 5 Compact Monitor	M-NESTPR-00-01	
DO06 P/N 2223	Datex-Ohmeda System 5 Compact Monitor	M-PRESTN..01	
DS01 P/N 1325	Corometrics Maternal / Fetal Monitor Model 118 Maternal / Fetal Monitor Model 120		
	Critikon Dinamap MPS Select Monitor		
	Datascope Expert Monitor with Nellcor SpO <sub>2</sub>	HS-500	
	Datascope Accutor Monitor with Nellcor SpO <sub>2</sub> Passport Monitor with Nellcor SpO <sub>2</sub>		
	North American Draeger Evita 2 Dura, Evita 4 Anesthesia Machine Julian Anesthesia Machine Narcomed 4 Anesthesia Machine Vitalert 1000, 3100 Monitor		
	Edentec / Nellcor Puritan Bennett Edentrace II Plus (Model 3711)		
	Medical Data Electronics Escort Prism Monitor	20001	
	Medical Data Electronics Escort / Escort II Monitor		
	Welch Allyn Protocol Propaq CS Monitor Propaq Encore Monitor Propaq Series 100 Monitor		
	Witt Biomedical Frontend Patient Monitor		
DS02 P/N 1539	Datascope Passport Monitor with Datascope SpO <sub>2</sub>		
	Datascope Accutor Monitor with Datascope SpO <sub>2</sub>		
DS03 P/N 1365	Datascope Expert Fukuda Denshi Dynascope DS-5300W		

**Note:** Do not interface the Radical with any monitor or module not specifically listed in this chart.

11683/2694S-0809


**\* Not for Sale in the United States**

SatShare™ Model	Validated SatShare™ Monitor Interfaces	Model / Plug-In Parameter Modules	Host Instrument Connector
HP03 P/N 1357	Philips / Agilent / HP Viridia CMS Model / Merlin CMS Model M1165, M1166, M1167, M1175, M1176, M1177 Monitors M1092, M1094, M1095, M1097 Display Modules, using the M1046 Computer Module (Main Processing Module)	M1020A	
	Philips / Agilent / HP Omnicare / Viridia Model 24/26 Series Monitor System M1204		
	Philips / Agilent / HP Model 78352 Monitor	M1020A	
	Philips / Agilent / HP Series 50XM Fetal Monitor		
MQ01 P/N 1321	GE / Marquette Solar 7000 / 8000 Series Monitor Tramscope 12 / 12C Monitor Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor Smart-Pac Transport Display Monitor	Tram 250SL Tram 450SL Tram 650SL Tram 850A&SL SpO <sub>2</sub> Solar Parameter Module	
	GE / Marquette Eagle 3000 Series Monitor		
	GE / Marquette Eagle 4000 Series Monitor		
	GE / Marquette Dash Series Monitor		
	North American Draeger Narkomed 6100 Anesthesia Systems with Integrated Patient Monitoring		
MS01 P/N 1789	Datex-Ohmeda S5	E-Masimo-00	
	Draeger Infinity Delta	SmartPOD with Masimo SET	
	GE/Marquette Solar 8000	Tram 451M Tram 851M Masimo SET single parameter module	
	Philips IntelliVue Monitor	MMS Module	
NK01 P/N 1367	Nihon Kohden Model BMS-8301A/MV-831RA, Model BMS-8502A, BSM2101, and BSM2102	AL-800PA	
NK02 P/N 1368	Nihon Kohden BSM-2301, BSM-2303, BSM-2351, BSM-2353, BSM-2391, BSM-4101, BSM-4103, BSM-4111, BSM-4113, BSM-5105, BSM-5132, BSM-5135, BSM-5192, BSM-9101, BSM-9800, and BSM-9502		
OH01 P/N 1366	Airshields Athena Monitor with Ohmeda SpO <sub>2</sub>		
	Ohmeda RGM Monitor (Respiratory Gas Machine)		
	GE / Marquette Solar 7000 / 8000 Series Monitor Tramscope 12 / 12C Monitor Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor Smart-Pac Transport Display Monitor	Tram 100A Tram 200A&SL Tram 400A Tram 500A Tram 600A&SL Tram 800A&SL	
OH02* P/N 1646*	Hellige SMU 611 Monitor*		
PM01 P/N 1840	Philips Medical M3/M4	MMS Module	
	Philips Intellivue Monitor	MMS Module	
SL01 P/N 1323	SpaceLabs Medical PC Express Monitor (90308) PC Scout Monitor (90309) PC Bedside / Central Monitor (90303, 90311, 90312) PC2 Bedside / Central Monitor (90305, 90315, 90316)	90465 90466 90467 90489	
	SpaceLabs Medical 600 Series Monitor		
SL02 P/N 1362	SpaceLabs Medical Ultraview 1030 / 1050 Portable Bedside Monitor (90367, 90369) Ultraview 1500 / 1600 / 1700 Bedside Monitor (90363, 90364, 90387) Ultraview UCW Bedside Monitor (90385)	90496	

**Note:** Do not interface the Radical with any monitor or module not specifically listed in this chart.

11683/2694S-0809

**\* Not for Sale in the United States**

SM01 P/N 1322	Siemens Medical Systems Sirecust 1281 Monitor SC 6002XL, SC7000, SC9000XL		
------------------	---	--	---

Monitor Manufacturer	Validated SatShare™ Monitor Interfaces	Validated SatShare™ Plug-In Parameter Modules	SatShare™ Model
Airshields	Athena Monitor with Nellcor SpO <sub>2</sub>		AT01
	Athena Monitor with Ohmeda SpO <sub>2</sub>		OH01
Corometrics	Maternal / Fetal Monitor Model 118		DS01
	Maternal / Fetal Monitor Model 120		DS01
Critikon	Dinamap MPS Select Monitor		DS01
Datascope	Accutor Monitor with Datascope SpO <sub>2</sub>		DS02
	Accutor Monitor with Nellcor SpO <sub>2</sub>		DS01
	Expert Monitor with Nellcor SpO <sub>2</sub>	HS-500	DS01
	Datascope Expert DS-5300W		DS03
	Passport Monitor with Datascope SpO <sub>2</sub>		DS02
	Passport Monitor with Nellcor SpO <sub>2</sub>		DS01
Datex-Ohmeda	AS/3 Patient Monitor	M-EST	DO03
	AS/3 Patient Monitor	M-ESTP	DO04
	AS/3 Patient Monitor	M-NESTPR-00-00	DO02
	AS/3 Patient Monitor	M-NESTR-00-00	DO02
	AS/3 Patient Monitor	M-NSAT	DO01
	CardioCap/5 Monitor		DO02
	Light Monitor		DO02
	S5	E-Masimo-00	MS01
	System 5 Compact Monitor	M-NESTPR-00-00 M-NESTR-00-00	DO02
	System 5 Compact Monitor	M-EST M-ESTP	DO04
	System 5 Compact Monitor	M-NESTPR-00-01	DO05
	System 5 Compact Monitor	M-NSAT	DO01
	System 5 Compact Monitor	M-PRESTN..01	DO06
Draeger	Infinity Delta	SmartPOD with Masimo SET	MS01
Dynascope	Dynascope DS-5300W		DS03
Edentec / Nellcor Puritan Bennett	Edentrace II Plus (Model 3711)		DS01
GE / Marquette	Dash Series Monitor		MQ01
	Eagle 3000 / Eagle 4000 Series Monitor		MQ01
	Solar 7000 / Solar 8000 Series Monitor Tramscope 12 / 12C Monitor Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor Smart-Pac Transport Display Monitor	Tram 250SL Tram 450SL Tram 650SL Tram 850A&SL SpO <sub>2</sub> Solar Parameter Module	MQ01
	Solar 7000 / Solar 8000 Series Monitor Tramscope 12 / 12C Monitor Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor Smart-Pac Transport Display Monitor	Tram 100A Tram 200A&SL Tram 400A Tram 500A Tram 600A&SL Tram 800A&SL	OH01
	Hellige Medical Systems SMU 611 Monitor*		OH02*
	Solar 8000	Tram 451M Tram 851M  Masimo SET single parameter module	MS01
Philips / Agilent / HP	Model 78352 Monitor	M1020A	HP03
	Omnicare / Viridia Model 24/26 Series Monitor System M1204		HP03
	Philips IntelliVue Monitor	MMS Module	MS01
	Philips IntelliVue Monitor	MMS Module	PM01
	Philips Medical M3/M4	MMS Module	PM01
	Series 50XM Fetal Monitor		HP03
	Viridia CMS Monitors	M1020A	HP03
Medical Data Electronics	Escort / Escort II Monitor		DS01
	Escort Prism Monitor	20001	DS01
Nihon Kohden	Model BMS-8301A/MV-831RA	AL-800PA	NK01
	Model BMS-8502A	AL-800PA	NK01
	Model BSM-2301A and BSM-5105K		NK02
North American Draeger	Evita 2 Dura, Evita 4 Anesthesia Machine		DS01

**Note:** Do not interface the Radical with any monitor or module not specifically listed in this chart.

11683/2694S-0809

**\* Not for Sale in the United States**

# SatShare™-Interface Cable

## Validation Chart

Monitor Manufacturer	Validated SatShare™ Monitor Interfaces	Validated SatShare™ Plug-In Parameter Modules	SatShare™ Model
	Julian Anesthesia Machine		DS01
	Narcomed 4 Anesthesia Machine		DS01
	Narkomed 6100 Anesthesia Systems with Integrated Patient Monitoring		MQ01
	<b>Vitalert 1000, Vitalert 3100 Monitor</b>		<b>DS01</b>
Ohmeda	RGM Monitor (Respiratory Gas Machine)		OH01
Siemens Medical Systems	Sirecust 1281 Monitor		SM01
	SC 6002XL, SC7000, SC9000XL		SM01
Spacelabs Medical	PC Express Monitor (90308) PC Scout Monitor (90309) PC Bedside / Central Monitor (90303, 90311, 90312) PC2 Bedside / Central Monitor (90305, 90315, 90316)	90465 90466 90467 90489	SL01
	600 Series Monitor		SL01
	Ultraview 1030 / 1050 Portable Bedside Monitor (90367, 90369) Ultraview 1500 / 1600 / 1700 Bedside Monitor (90363, 90364, 90387) Ultraview UCW Bedside Monitor (90385)	90496	SL02
Welch Allyn Protocol	Propaq CS Monitor		DS01
	Propaq Encore Monitor		DS01
	Propaq Series 100 Monitor		DS01
Witt Biomedical	Frontend Patient Monitor		DS01

**Note:** Do not interface the Radical with any monitor or module not specifically listed in this chart.

11683/2694S-0809

**\* Not for Sale in the United States**