

Noninvasive Blood Pressure and Temperature Accessories

For Use with Root® Patient Monitoring and Connectivity Platform and Rad-97™ Pulse CO-Oximeter®



Root with Noninvasive Blood Pressure and Temperature

Rad-97 with Noninvasive Blood Pressure






Masimo Blood Pressure Cuffs

- > Disposable and reusable blood pressure cuffs, for a variety of patient types, designed for reliability and patient comfort
- > Interior and exterior indicators facilitate proper cuff size selection and fit
- > 100% latex-free materials
- > Meets AHA, AAMI-SP10, and MDD guidelines

Disposable Blood Pressure Cuffs




- > Single-patient-use
- > Durable, single-piece design allows for multiple inflations
- > For use with Root NIBPT and Rad-97 NIBP

Neonatal:

Part #	Description	Range	Color	Cuff
25242	Neonate Size 1 Cuff (20/box) Male Bayonet, for use with Female Bayonet NIBP Hose	3-6 cm	Orange	
25243	Neonate Size 2 Cuff (20/box) Male Bayonet, for use with Female Bayonet NIBP Hose	4-8 cm	Turquoise	
25244	Neonate Size 3 Cuff (20/box) Male Bayonet, for use with Female Bayonet NIBP Hose	6-11 cm	Green	
25245	Neonate Size 4 Cuff (20/box) Male Bayonet, for use with Female Bayonet NIBP Hose	7-13 cm	Navy Blue	
25246	Neonate Size 5 Cuff (20/box) Male Bayonet, for use with Female Bayonet NIBP Hose	8-15 cm	Burgundy	

Blood Pressure Hoses:

- > For NIBP measurement, Root NIBPT and Rad-97 NIBP come with a blood pressure hose (choice of male quick connect or female bayonet, and adapter if needed). Blood pressure cuffs can be ordered separately.

Part #	Description	Cuff Connector	Hose
25247	NIBP Patient Hose, 3 Meters, Male Quick Connect	Female Quick Connect	
25302	NIBP Patient Hose, 3 Meters, Female Bayonet	Male Bayonet	
26426	NIBP Patient Hose Adapter, 10 cm, Male Bayonet to Female Slip Luer (For use with GE CRITIKON Neonatal BP cuffs)	Male Slip Luer	

Adult & Pediatric:

Part #	Description	Range	Color	Cuff
25296	Disposable Infant Cuff (20/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	8-13 cm	Orange	
25297	Disposable Child Cuff (20/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	12-19 cm	Green	
25298	Disposable Small Adult Cuff (20/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	17-25 cm	Turquoise	
25299	Disposable Adult Cuff (20/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	23-33 cm	Navy Blue	
25300	Disposable Large Adult Cuff (20/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	31-40 cm	Burgundy	
25301	Disposable Thigh Cuff (20/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	38-50 cm	Brown	


Reusable Blood Pressure Cuffs

- > Simple to clean and disinfect without compromising longevity
- > Unique tapered tail design allows the cuff to be easily wrapped and applied to the patient
- > Durable, single-piece design allows for use on multiple patients
- > 100% PVC-free and DEHP-free
- > For use with Root NIBPT and Rad-97 NIBP

Adult & Pediatric:

Part #	Description	Range	Color	Cuff
25303	Reusable Infant Cuff (5/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	8-13 cm	Orange	
25260	Reusable Child Cuff, Masimo (5/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	12-19 cm	Green	
25261	Reusable Small Adult Cuff, Masimo (5/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	17-25 cm	Turquoise	
25262	Reusable Adult Cuff, Masimo (5/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	23-33 cm	Navy Blue	
25263	Reusable Large Adult Cuff, Masimo (5/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	31-40 cm	Burgundy	
25304	Reusable Thigh Cuff (5/box) Female Quick Connect, for use with Male Quick Connect NIBP Hose	38-50 cm	Brown	

Blood Pressure Hose:

Part #	Description	Cuff Connector	Hose
25247	NIBP Patient Hose, 3 Meters, Male Quick Connect	Female Quick Connect	

Compatibility and Warranty Information

Compatibility

In addition to the NIBP cuffs and hoses offered by Masimo, Root NIBPT and Rad-97 NIBP are also compatible with GE CRITIKON SOFT-CUF Neonatal blood pressure cuffs when utilizing the patient hose adapter.

- > GE CRITIKON Part Numbers: 2121, 2122, 2124, 2125, 2126. (not available for order from Masimo)

The products listed in this catalog and GE CRITIKON SOFT-CUF Neonatal cuffs have been validated for use with the NIBP technology offered in Root NIBPT and Rad-97 NIBP. Additional blood pressure cuffs and hoses provided by other manufacturers have not been validated with the technology, and therefore the performance and compatibility cannot be guaranteed.

Warranty

SunTech Directions for Use (DFU) state:

Reusable Cuffs: SunTech Medical warrants the blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. This limited warranty covers the no charge replacement of the cuff under normal wear and tear conditions when returned to the attention of Service Department at the address below depending on location.

Disposable Cuffs: SunTech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. The warranty shall only cover out-of-box failures, as this is a Single Use product. It is not intended for repeated use.

Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Service Department at one of the addresses below.

SunTech Medical, Inc.

507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200, USA
TEL: 1.800.421.8626 or 1.919.654.2300
FAX: 1.919.654.2301

SunTech Medical, Ltd. Europe



Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham,
Oxfordshire OX29 4TS
England
TEL: 44.1 865.884 234
FAX: 44.1 865.884 235

Temperature Probes

> For use with Root NIBPT

- For temperature measurement, Root NIBPT comes with a thermometer probe well kit (choice of oral/axillary or rectal). Temperature probe covers can be ordered separately.

Adult & Pediatric:

Part #	Description	Temperature Accessories
25251	Oral Temperature Probe and Well Kit, Welch Allyn	
25252	Rectal Temperature Probe and Well Kit, Welch Allyn	
26289	Temperature Probe Covers, Welch Allyn (7500/case)	

Additional Accessories

Part #	Description	Image
24205	Root Roll Stand	
19118	Barcode Reader, Black	
19274	Root Roll Stand Scanner Bracket Accessory	
38532	Root Scanner Wall Mount Accessory	
24209	Roll Stand Cavi Wipes Bracket Accessory (Compatible with roll stands for Root and Rad-97)	
38895	Rad-97 Roll Stand	
19356	Locking Power Cord, US, 10A, 125V, RoHS	
18313	USB-Null Modern Cable (2.5m)	

*The use of the trademark PATIENT SAFETYNET is under license from University Health System Consortium.

Rad-97 NIBP is not licensed for sale in Canada.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

